

PREPARATION / IRON DECONTAMINATION

IRON

NL-TDS-IRON-EN | REV 1.0 | DOCUMENT DATE 26 SEP 2026

ENGLISH EDITION / PRODUCT TECHNICAL INFORMATION

An iron-oxide remover with slow-drying, low-odour behaviour. A purple reaction indicates contact with ferrous contamination.

Intended surfaces: Painted surfaces, wheels, chrome, steel and aluminium, subject to product instructions.

01 / TECHNICAL DATA

Pack size	473 mL / 16 fl oz
Typical reaction time	Approx. 3 min
Reaction indicator	Purple colour
Application method	Spray; rinse thoroughly
Surface condition	Cool; not dripping wet
pH / dilution	Not specified

Technical values refer to the product information [1]. Where stated, curing and working times are typical guidance and depend on conditions.

04 / STORAGE & SHELF LIFE

Store at 0-25°C in the tightly closed original container. Shelf life: at least 12 months; use within 180 days after opening.

05 / APPLICATION NOTES & SAFE USE

The purple reaction is a visual application indicator. Follow the recommended dwell time and rinse all residues thoroughly.

For technical parameters not specified in this document and substrate-specific advice, contact NANOLAB technical support.

Read the current product-specific SDS and label before use. Follow their ventilation, PPE, handling and disposal instructions. This sheet does not replace the SDS.

[1] Product information: www.nanolab.ltd/iron/

Reference information accessed 26 Sep 2026. Use together with the current product label and SDS.

02 / APPLICATION

01 Wash; remove excess water. Work on a cool surface with ventilation, gloves and eye protection.

02 Spray onto contamination; allow about 3 minutes for reaction.

03 Gently agitate if needed, rinse thoroughly, inspect and repeat if necessary. Dry with microfiber.

03 / CARE & MAINTENANCE

Use as a decontamination stage when required.